

OvationMedical<sup>®</sup>  
**The Game Changer**<sup>™</sup> *Gen2*



Premium Universal OA Knee Brace

**INSTRUCTIONS FOR USE**

Ovation Medical®  
**The Game Changer™ Gen2**

Simple Adjustment Hinge



Buckle / Strap

#4

Buckle / Strap

#3

Condyle Pad

Buckle / Strap

#1

Buckle / Strap

#2

Adjustment  
Key



For assistance please contact our Provider/Patient OA support line  
Call toll free **1.844.357.4931** M-F 7am-3pm PST  
[ovationmedOA.com](http://ovationmedOA.com)



## INITIAL FITTING PROCESS



1. Release all buckles by pushing the buckle release buttons. Detach clasps from frame



2. Position brace by sitting on stable surface with knee bent at approx. 45 degrees. Align the center of the hinge to just above the patella



3. Reattach the straps in sequence (1-4). Attach all buckles in sequence. If straps need further adjustment see page 5




4. For final brace adjustment, verify proper hinge alignment (Step 2) then adjust strap tabs

For assistance please contact our Provider/Patient OA support line  
Call toll free **1.844.357.4931** M-F 7am-3pm PST  
**ovationmedOA.com**

## MEDIAL / LATERAL OFF-LOADING

APPLY BRACE ON LATERAL (OUTSIDE) OF KNEE

### MEDIAL OFF-LOADING




MEDIAL  
COMPARTMENT  
OFF-LOADING  
(PAIN RELIEF)



**TURN KEY  
COUNTER  
CLOCKWISE  
(IN)**

### LATERAL OFF-LOADING



LATERAL  
COMPARTMENT  
OFF-LOADING  
(PAIN RELIEF)



**TURN KEY  
CLOCKWISE  
(OUT)**

For assistance please contact our Provider/Patient OA support line  
Call toll free **1.844.357.4931** M-F 7am-3pm PST  
[ovationmedOA.com](http://ovationmedOA.com)



## BRACE COMFORT ASSESSMENT



Once the brace is applied, stand and walk a short distance to ensure comfort and proper correction. If there is any discomfort, adjust the straps as needed. If the patient experiences pain in the knee joint or leg, turn the hinge adjustment back slightly. If more relief is required turn the hinge adjustment no more than 1/4 of a turn; as directed by your provider.




## SUSPENSION WRAP



Apply Suspension Wrap to lower portion of brace (as shown) for improved fit. Do not over tighten.

For assistance please contact our Provider/Patient OA support line  
Call toll free **1.844.357.4931** M-F 7am-3pm PST  
[ovationmedOA.com](http://ovationmedOA.com)

# STRAP ASSEMBLY / ADJUSTMENT



1. Determine desired strap length



2. Remove strap pads



3. Remove strap ends



4. Detach velcro tab from strap



5. Adjust strap length, trim strap as needed



6. Reattach strap ends and reapply strap pads

For assistance please contact our Provider/Patient OA support line

Call toll free **1.844.357.4931** M-F 7am-3pm PST

[ovationmedOA.com](http://ovationmedOA.com)



Ovation Medical®  
**The Game Changer™ Gen2**  
**Wear Time**

Start with 30 minutes on day one and gradually increase brace wear time to allow your body to adjust to the brace.

Patients can expect clinically relevant pain relief and improvement of knee function over time.

Our slim fitting brace design can be used under pants or shorts.



### PROVIDER/PATIENT OA SUPPORT LINE

Ovation Medical is fully committed to support both for providers and patients to get the best possible results with your OA knee Brace. Healthcare providers AND patients should feel free to call the Provider/Patient OA support line for any reason. This may include fitting assistance, instructions for wear, wash instructions, or any question or support you require. Our Provider/Patient OA support line is available Monday through Friday, 7am - 3pm at **1.844.357.4931**



For assistance please contact our Provider/Patient OA support line  
Call toll free **1.844.357.4931** M-F 7am-3pm PST  
**ovationmedOA.com**

## FAQS

Q: What if I still cannot get the brace to fit properly?

A: Contact your provider or feel free to contact our Provider/Patient OA support line at 1.844.357.4931

Q: Can I get my brace wet?

A: Yes, the hinges will not rust. Be sure to rinse and air dry the straps and liners.

Q: If a strap needs to be replaced, what do I do?

A: Contact your provider for a replacement strap kit.

## CARING FOR YOUR BRACE

1. To clean the liners, straps, pads and suspension wrap:
  - A. Remove from brace.
  - B. Hand wash with cold water and MILD detergent.
  - C. Rinse well with cold water.
  - D. Let it air dry; Clean & dry liners will result in less skin irritations.
2. To clean the brace frame and hinge; wipe with a damp cloth.
3. Avoid brace contact with salt water, sand, dirt and chlorinated water.
4. Use a dry silicone spray if a hinge lubricant is needed.

**For warranty information:**

Please visit our website at [ovationmedOA.com](http://ovationmedOA.com)

The Game Changer Copyright 2017