EN PREMIUM UNIVERSAL OA KNEE BRACE.
THE GAME CHANGER®





DE UNIVERSAL-PREMIUM-OA-KNIESCHIENE.
THE GAME CHANGER®

ES RODILLERA PREMIUM PARA ARTROSIS,
TALLA ÚNICA. THE GAME CHANGER®

FR ATTACHE POUR GENOU OA UNIVERSEL
PREMIUM. THE GAME CHANGER®

IT GINOCCHIERA TAGLIA UNICA PREMIUM OA.
THE GAME CHANGER®

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ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The Premium Universal OA Knee Brace is intended for single compartment off-loading of an osteoarthritic knee. The Premium Universal OA Knee Brace provides support for patients suffering from Osteoarthritis of the knee.

The Premium Universal OA Knee Brace is intended for short term single patient use and not intended to undergo regular service cycles.



This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

CONTRAINDICATIONS: None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

WARNINGS AND PRECAUTIONS: Do not over tighten, adjust as necessary. If swelling increases and pain, numbness, skin irritation occurs, discontinue use. Contact medical provider immediately. This product is only a support device, not to prevent injuries or prevent recurrence of injuries. This product is to be used under the supervision of a medical professional.

REPORTING PROBLEMS: The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

PRODUCT DESCRIPTION:

APPLICATION INSTRUCTIONS:

- Release all buckles by pushing the buckle release buttons. Detach clasps from frame.
- Position brace by sitting on stable surface with knee bent at approx. 45 degrees. Align the centre of the hinge to just above the patella.
- Reattach the straps in sequence (1-4). Attach all buckles in sequence. If straps need further adjustment see page 5.
- For final brace adjustment, verify proper hinge alignment (step 2) then adjust strap tabs.

STRAP ASSEMBLY/ADJUSTMENT:

- Determine desired strap length.
- Remove strap pads.
- Remove strap ends.
- Detach velcro tab from strap.
- Adjust strap length, trim strap as needed.
- Reattach strap ends and reapply strap pads.

BRACE COMFORT ASSESSMENT: In the event of discomfort adjust straps and hinge and/or contact your provider. Once the brace is applied, stand and walk a short distance to ensure comfort and proper function. If there is any discomfort, adjust the straps as needed. If the patient experiences pain in the knee joint or leg, turn the hinge adjustment back slightly. If more relief is required turn the hinge adjustment no more than 1/4 of a turn; as directed by your provider.

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Si desea información sobre la garantía, visite nuestro sitio web en la dirección ovationmedOA.com. The Game Changer Copyright 2017.
Pour toute information sur la garantie, veuillez consulter notre site Web à l'adresse ovationmedOA.com. The Game Changer Copyright 2017.
Per informazioni sulla garanzia, visitate il nostro sito web all'indirizzo www.ovationmedOA.com. The Game Changer Copyright 2017.