EN COMPACT PRO POST-OP KNEE

Quick Release Buckle, Cool Wrap with Drop Lock and Extending Struts

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







**1.****2.****3.****4.****5.****6.****7.****8.****ENGLISH**

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The Compact Pro Post-Op Knee is intended for recovery of post-operative of a knee and it is indicated for ACL/PCL/MCL/LCL, osteochondral, patella tendon and meniscal repairs, tibial plateau and condylar fractures, sprains and strains of the knee as well as HTOs (High Tibial Osteotomy) and other injuries.

The Compact Pro Post-Op Knee is intended for short term single patient use and not intended to undergo regular service cycles.

This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

CONTRAINDICATIONS: None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

WARNINGS AND PRECAUTIONS: If pain and/or numbness occurs, contact your medical professional immediately. Always consult with your physician or therapist before making changes to the settings. This product is to be used under the supervision of a medical professional. This device is not intended for unsupervised public use. This product is designed for single patient use only.

REPORTING PROBLEMS: The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

APPLICATION INSTRUCTIONS:

1. Open brace as shown.
2. Position brace in a flat position as shown. Set leg over brace centered approximately between both struts.
3. Leg should be positioned approximately as shown centered between the struts with straps positioned to the rear of the leg. Verify hinge and buckle orientation is as shown.
4. Each side of the brace has 3 sliders that can be independently adjusted to accommodate different leg lengths. To adjust push slider button and move to desired position. Ensure that the button 'snaps' and locks into new position for secure placement and that both sides of brace are set at the same positions.
5. Wrap brace around both sides of leg. Struts should be positioned as shown in the side view, with hinge slightly above center of knee. Secure straps in the following sequence: 3, 2, 4, 1 by inserting quick release buckle and tightening with free end. Excess strap ends can be trimmed off.
6. ROM (Range of Motion) dials are unlocked by pushing button toward center of hinge. Push button and rotate ROM dials to desired flexion/extension position and release button to lock. After adjustment make sure both left and right sides of brace are set at the same positions.
7. The brace may be locked into one of 4 positions (0, 15, 30 and 45 degrees). Push Drop Lock lever DOWN to lock and UP to unlock. Always make sure both left and right locks are set at the same positions.
8. If desired the ROM Adjusters can be locked in position with included zip ties. To lock insert zip tie as shown making sure it is behind button tab with blue printed arrow. Wrap around ROM adjuster and secure snugly. Trim off excess after tightening.

CARE AND MAINTENANCE: Hand wash in cold water using mild soap, do not bleach, line dry only, avoid direct heat or sun. Store in a cool dry place.

DEVICE DISPOSAL: Ovation Medical® products should be properly disposed of by a qualified recycler when the post-operative use period is completed or when the device is no longer required by the patient.

SYMBOLGY:

- a. Consult instructions for use and accompanying documents.
- b. Manufactured by. Note: can be used in combination with date of manufacture.
- c. Date of manufacture accompanied by a date.
- d. Single patient multiple use.
- e. Device sold in non-sterile condition.
- f. Not made with Natural Rubber Latex
- g. Peel here.
- h. Lot number.
- i. Reference number or product number.
- j. Quantity contained in package.
- k. UDI Unique Device Identifier Number.
- l. Product intended to be used as a medical device.
- m. Authorized representative (REP) in the European Community (EC).
- n. Conform to the European Medical Device Regulation (EU) 2017/745 for class I medical devices

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