

Dorsal Night Splint

APPLICATION INSTRUCTIONS:

- 1. Undo toe and calf straps.
- 2. Step foot into brace.
- 3. Secure toe strap by affixing the contact closures.
- 4. Apply calf strap by wrapping behind the leg and affix contact closures.

Washing Instructions:

Hand wash in cold water using mild soap. Do not bleach. Line - dry only. Avoid direct heat or sun.

NOTE: Do not over tighten, adjust as necessary. **CAUTION:** If swelling increases and pain, numbness, skin irritaion occurs, discontinue use. Contact medical provider immediately. This product is only a support device, not to prevent injuries or prevent reoccurance of injuries.



CE

Rev. 1