



Neoprene Knee Support with Stabilized Patella

APPLICATION INSTRUCTIONS:

1. Unfasten straps.
2. Place foot into the widest opening of the brace.
3. Pull the brace up the leg until the opening is directly over the knee cap.
4. Loop straps through the rings and pullback over the top of the lower thigh and upper calf fastening the hook and loop closure.
5. Adjust straps for comfort and fit (The brace should be snug but not uncomfortably tight).
6. After adjustment the straps may be cut (Be sure and leave enough strap for readjustment - approx. 5").

NOTE: Do not over tighten, adjust as necessary.

CAUTION: If swelling increases and pain, numbness, skin irritation occurs, discontinue use. Contact medical provider immediately.

This product is only a support device, not to prevent injuries or prevent reoccurrence of injuries.



OvationMedical® EC REP

5155 Clareton Dr.
Agoura Hills, CA. 91301
USA
800 403-6466
www.ovationmed.com

MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany