- EN** NU-FORM® UNIVERSAL BACK BRACE 0627
- DE** NU-FORM® UNIVERSAL-RÜCKENSTÜTZE 0627
- ES** ORTESIS DORSAL DE TALLA UNIVERSAL NU-FORM® 0627
- FR** ATTELLE POUR DOS UNIVERSELLE NU-FORM® 0627
- IT** FASCIA POSTURALE PER LA SCHIENA TAGLIA UNICA 0627
- PL** UNIWERSALNA ORTEZA KRĘGOSŁUPA 0627 NU-FORM®

**NU-FORM™**

**ENGLISH**

**BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.**

**INTENDED USE/INDICATIONS:** The Nu-Form® Universal Back Brace is intended to provide support for the lumbar region and it is indicated for lower back pain, acute lumbar sprains and strains, post-surgical fusion, laminectomy, discectomy, congenital or traumatic instability, disc herniation, spinal stenosis, postural backache, fractures, and degenerative joint/disc disease.

The Nu-Form® Universal Back Brace is intended for short term single patient use and not intended to undergo regular service cycles.


This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

**CONTRAINDICATIONS:** None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

**WARNINGS AND PRECAUTIONS:** Do not over tighten, adjust as necessary. If swelling increases and pain, numbness, skin irritation occurs, discontinue use. Contact medical provider immediately. This product is only a support device, not to prevent injuries or prevent recurrence of injuries. This product is to be used under the supervision of a medical professional.

**REPORTING PROBLEMS:** The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

**PACKAGE CONTENTS:**

**APPLICATION INSTRUCTIONS:**

1. Position the brace so that it is centered on the back over LS area.
2. Pull the belt panel ends forward around waist, placing right belt panel on top of left belt. Engage the hook and loop.
3. Pull both sides of pull rings forward and around for proper fit (micro fit).
4. Secure pull rings by pressing firmly onto front of the brace. Adjust as needed. Security strap may be used for additional control.

**ADJUSTMENT INSTRUCTIONS:**


Fig. 2: Adjust sizing by moving panels through D-ring to the desired size. Size markings are clearly labeled on the inside of each panel. Align ends of panel to desired size indicators then firmly secure each panel by attaching the hook to loop area.


Fig. 3: Attach Anterior Panel (#2) to the inside of the left belt panel.

**EXTENSION PANEL APPLICATION:**


Fig. 4: Extension Panel application. If patient requires a larger size than universal, an Extension Panel is available to add 18 inches in length (70 inches in total girth). The Extension Panel attaches to the outside of the left belt panel. The Velcro flap of the Extension Panel should attach to the inside of the belt panel. Replace Anterior Panel (#2) as shown in Fig. 3.

**CARE AND MAINTENANCE:** Hand wash in cold water using mild soap, do not bleach, line dry only, avoid direct heat or sun. Store in a cool dry place.

**DEVICE DISPOSAL:** Ovation Medical® products should be properly disposed of by a qualified recycler when the post-operative use period is completed or when the device is no longer required by the patient.

**SYMBOLGY:**

- a. Consult instructions for use and accompanying documents.
- b. Manufactured by. Note: can be used in combination with date of manufacture.
- c. Date of manufacture accompanied by a date.
- d. Single patient multiple use.
- e. Device sold in non-sterile condition.
- f. Not made with Natural Rubber Latex
- g. Peel here.
- h. Lot number.
- i. Reference number or product number.
- j. Quantity contained in package.
- k. UDI Unique Device Identifier Number.
- l. Product intended to be used as a medical device.
- m. Authorized representative (REP) in the European Community (EC).
- n. Conform to the European Medical Device Regulation (EU) 2017/745 for class I medical devices

**THIS PRODUCT MAY BE SUBSTANTIALLY MODIFIED BY AN INDIVIDUAL WITH EXPERTISE**  
**For more information, please contact customer care (U.S only) (844) 357-4931 or reach out to your Ovation Medical® representative.**

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