EN	NU-FORM® ANKLE BRACE
DE	NU-FORM® KNÖCHELSTÜTZE
ES	TOBILLERA NU-FORM®
FR	ATTELLE DE CHEVILLE NU-FORM®
IT	CAVIGLIERA CON CINGHIE NUFORM®
PL	ORTEZA STAWU SKOKOWEGO NU-FORM®

NU-FORM™

ENGLISH

BEFORE USING THE DEVICE, PLEASE READ THE FOLLOWING INSTRUCTIONS COMPLETELY AND CAREFULLY. CORRECT APPLICATION IS VITAL TO THE PROPER FUNCTIONING OF THE DEVICE.

INTENDED USE/INDICATIONS: The Nu-Form® Ankle Brace intended use is as a support for lower leg, ankle, and toe region as well as prophylactic support of foot and ankle during full activity. The product is indicated for plantar fasciitis, Achilles tendonitis, drop foot, post-static pain as well as inversion or eversion ankle sprains; medial or lateral ankle instability and subtalar joint instability.

The Nu-Form® Ankle Brace is intended for short term single patient use and not intended to undergo regular service cycles.

This product was designed to withstand normal use conditions for orthopedic support, post injury rehab or prevention of further injury. Like all products that are worn by individuals, the extended life of such products is determined by the level of activity and type of activity of the individual.

CONTRAINDICATIONS: None known. In every case, the final decision whether to use this Ovation Medical® device rests with the attending physician/therapist. He/she should assess prior to use of the device, if it is suitable for a given application, taking into account the condition of the patient.

WARNINGS AND PRECAUTIONS: Do not over tighten, adjust as necessary. If swelling increases and pain, numbness or skin irritation occurs, discontinue use and contact medical provider immediately. This product is only a support device, not to prevent injuries or prevent reoccurrence of injuries. This product is to be used under the supervision of a medical professional.

REPORTING PROBLEMS: The user should report any suspected serious incident related to the use of this Ovation Medical® device by informing the manufacturer and the competent authority of the state or the country in which the incident occurred.

APPLICATION INSTRUCTIONS:

1. Unfasten the Velcro straps and untie the laces.
2. Loosen the laces enough to enable the foot to slide easily into the brace.
3. Position heel into heel opening.
4. Tighten the laces snugly. Tie the laces.
5. Apply the figure-8 straps. Pull snugly.
6. Apply each strap to the Velcro tabs on each side of the ankle brace.

CARE AND MAINTENANCE: Hand wash in cold water using mild soap, do not bleach, line dry only, avoid direct heat or sun. Store in a cool dry place.

DEVICE DISPOSAL: Ovation Medical® products should be properly disposed of by a qualified recycler when the post-operative use period is completed or when the device is not longer required by the patient.

SYMBOLOLOGY:

- a. Consult instructions for use and accompanying documents.
- b. Manufactured by. Note: can be used in combination with date of manufacture.
- c. Date of manufacture accompanied by a date.
- d. Single patient multiple use.
- e. Device sold in non-sterile condition.
- f. Not made with Natural Rubber Latex
- g. Peel here.
- h. Lot number.
- i. Reference number or product number.
- j. Quantity contained in package.
- k. UDI Unique Device Identifier Number.
- l. Product intended to be used as a medical device.
- m. Authorized representative (REP) in the European Community (EC).
- n. Conform to the European Medical Device Regulation (EU) 2017/745 for class I medical devices

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